

Quazar IPL280

Instruction Material for Quazar IPL280
Precision O.E.M. IPL and Driver Unit

Quazar IPL and Driver

Replacement Parts

IPL Supplies

IPL Eyewear

Eyewear # D-213-4600 Each \$199.95

IPL280 Lamp Replacement

IPL # D-213-6700 Each \$299.95

Carbon Dye 200ml

Carbon Dye # P-213-2500 Each \$99.95

Prices are subject to change without notification. To order on-line, go to
<http://www.centre-biotechnique-avance.com>

For technical assistance beyond what this manual provides, please e-mail

admin1@centre-biotechnique-avance.com

Please allow 24 hours for processing.

Quazar IPL produces laser radiation which can be harmful to the eyes. Always wear protective eyewear while operating this equipment. Intense Pulsed Light has the capability to burn the skin if the technician does not closely observe the patient's reaction to the procedure. IPL treatments result in full destruction of the hair follicle and is **irreversible**. Always plan ahead before undertaking detail work such as eyebrow shaping or hairline contouring. Patch test a small area (no larger than 1X1 inch square) before full application. Allow 24 hours to determine the patient's reaction before applying full treatment.

Handling Precautions: The IPL flash tube can be damaged beyond repair if the handpiece is dropped or bumped sharply. Handle with care.

Quick Setup Guide


Read this guidebook first to set up your
equipment for use.



Keep this manual in a convenient place for quick and easy
reference at all times.

The product names in this guidebook are trademarks or registered trade marks of each specific manufacturer. In the interest of providing superior equipment, Quazar Industries reserves the right to modify or amend equipment specifications without notice or obligation.

IPL280 Instrument



The IPL280 comes with a 54fcm2 adjustable pulse frequency (fixed pulse duration) instrument for superior operator control. Dropping or bumping the instrument may result in irreversible damage to the internal components and would not be covered by warranty. Activating the IPL while the glass condenser is in contact with an object (such as the carbon test sheet) may damage the optics. When testing against the carbon sheet always hold the condenser at least 5mm from the surface.

Trouble Shooting


Should you encounter technical problems with your Quazar IPL system, refer to the following guide for potential problems and their solutions.

- Unit is plugged into the wall, all accessories are correctly inserted into the unit but no laser output is being registered.
- ++Check all connections. Plug and unplug each one being sure all contacts are sound.
- ++Check all cords. Due to continual bending and fatigue, wires may fray or break resulting in full loss of power.
- ++Check fuse: The Quazar unit has a fast-acting fuse mounted on the top panel of your unit to protect the delicate flash tube from voltage spikes, on-line power surges and electro-static discharge (ESD damage). Replace with 15 amp fast acting type only. Failure to comply with these specifications may result in serious damage to your IPL and will void all warranties.
- Unit hums or makes noises.
- ++Unit needs servicing.
- IPL output is weak.
- ++emitter optic is blocked with carbonized debris (burnt hair etc.). Clean thoroughly with a cotton tipped applicator and alcohol. If the performance of the IPL does not improve after cleaning the head your unit needs servicing by a qualified agent. Contact technical support for assistance.
- No output from the IPL is registered after all trouble-shooting suggestions listed above have been checked.
- ++Probable flash tube failure. Unit needs servicing.




Digital Temperature Meter

The meter may need re-calibration by adjusting the on-board potentiometer. Two resistors, Ra and Rb may be used in order to alter the full scale reading (f.s.R.) of the meter.



Specification	Min	Typ.	Max.	Unit	Reading
Accuracy		0.05	0.1	%(+1 count)	50
Linearity			+1	count	40
Sample rate			3	Samples/sec	30
Op. temp	0	50		°C	20
Temp. Stbl.		150		Ppm/°C	10

IPL280 Control Locations/Feature Descriptions



- A. **Key Lockout:** This feature is required by law on all high-power laser and light devices. The first step in the correct sequence to power-up your laser is to turn this clockwise using the special key included in your kit. The light (center LED lamp) will be off when the power system is off and red when the power system is on.
- B. **System Status LED:** This light will indicate system status (neutral or live). Blue indicates your IPL is live and ready for use.
- C. **Pulse Control:** This allows the operator to set the number of laser pulses per second.
- D. **IPL Jacks:** An eight-prong coaxial power jack for the diode laser hand piece and two blue banana plug receptacles.
- E. **Power Cord:** Rated for 60Hz, 120-240 V, 10 Amp with corresponding plug-style for country of destination.
- F. **IPL Instrument:** A 280 Watt high intensity pulsed light with thumb switch.
- G. **Eyewear:** This is an essential part of the treatment process. Direct or reflective IPL radiation can seriously injure the eye. Both the technician and the patient must use the protective eyewear while the IPL is enabled or activated. Eyewear is intended for **accidental** exposure only. Never stare directly into an IPL optic.
- H. **Carbon Dye:** This is an 'atomized' form of molecular carbon which easily penetrates deeply into the follicle shaft. The dye adds pigment which gives a receptor for the photon/heat exchange reaction. The carbon atoms will capture the IPL energy and convert it into heat for the rapid and efficient cauterization of tissue for the permanent destruction of the hair follicle organ.
- I. **High-Precision Tweezers:** Apparatus for the extraction of follicle prior to carbon dye application.
- J. **Temperature Readout:** Control feature for the regulation of the flash tube internal temperature. Do not exceed 50C.
- K. **Precision Optics:** BK7 Glass condenser optic for the efficient focus and penetration of light energy through the skin. May be changed to ruby or sapphire (see supplies).

Equipment Warranty

We warrant to the original purchaser the equipment manufactured by us to be free from defects in material and workmanship under normal use and service. Our obligation under this warranty shall be limited to the repair or exchange of any part or parts which may prove defective under normal use and service within 12 calendar months from the date of shipment and which our examination shall disclose to our satisfaction to be thus defective. When necessary, purchaser shall apply for a Return Materials Authorization and instructions on proper return procedures from their original sales associate. The laser diode (neod) requires special operating precautions which, if defied, may void warranty.

Warranty Extension Certification:

Customer Number	Authorization Number
Warranty Extension () years	Warranty Type: A B C D


IPL Startup Procedure

Put on your IPL eyewear before powering up your system. Wear the eye protection throughout the entire procedure.



Place the blue plugs into the IPL driver unit firmly. These jacks produce high voltage current so the system must be off during plug insertion (key lock in the 'lock off' position). As an added safety measure, do not touch the metal portion of the plug while inserting or removing from the driver (even if you assume the unit is off).

Turn the key lock switch from system status neutral (off) to system status blue (enabled). At this time your IPL hand piece is LIVE. Extreme care must be exercised while handling the IPL. This system produces high levels of potentially hazardous light which could injure unprotected eyes.



Testing the IPL

With the pulse setting on 6 (seven pulses per second), place the IPL condenser optic 5mm from the black test sheet included with your instruction material. Press the trigger switch then immediately release (just for testing purposes). The hand piece will deliver 7 pulses of light. Repeat the test but hold the trigger switch down for 10 seconds. This will deliver approximately 70 pulses which will ignite the test pad producing vapor and smoke.

It is recommended to use the IPL in 6 pulses per second (pps) at all times (setting 4) to avoid over treatment.

IPL Safety Considerations

IPL should not be used for actual 'IPL (i.e. laser) shaving'. The energy emission exit port (glass optic) is in direct contact with the skin during the treatment. The reaction that occurs when radiation strikes a hair is very intense heat. It could cause pits and defects in the optic as well as potentially burn the patient's skin which is in direct contact with the vaporizing hair.

IPL must be used on shaved (or waxed) skin. The follicle being targeted is about 1/8 to 1/4 inch below the skin (the follicle and papilla). If using carbon dye, remove any excess which may come in contact with the optics using ethyl alcohol.




Lamp Overheating

Keep the flash tube temperature at or below 50 degrees Centigrade for optimal performance. Pushing the IPL beyond the limits of performance (operating above 50C) may damage the tube, resulting in loss of intensity, performance and life expectancy of the IPL (which is not covered by your standard warranty).


IPL Optics Maintenance

Clean the condenser optics frequently. For light maintenance use acetone or alcohol on a cotton tip applicator. For more aggressive cleaning use an 'Exacto' blade to scrape any carbonized material from the glass.



Pulse Counter

Your IPL280 comes with a convenient pulse counter to track total photon emissions (for conversion to joules, see back page for information). To reset to zero press the button to the right of the display. The system will pulse 2 times per second on setting 2, 3 times on 3, 4 times on 4, 5 times on 5, and 7 times on 6.



Shipping Damage: The high wattage flash tube in your hand piece is sensitive to shock and may be damaged if the unit is dropped or mishandled. If your system arrives in non-working condition (the hand piece does not produce pulses of light energy when activated) it has likely been damaged during shipping. See 'Lamp Replacement Procedure' for instructions on removal and replacement of the high wattage tube.

Intense Pulsed Light Progressively Permanent Hair Removal

The IPL can be considered an excellent replacement for laser hair removal as the results and overall procedure are similar in many ways. It can be noted that IPL treatments are more efficient due to the substantial increase in overall energy being used (which equates to less overall treatment time).



Like laser, IPL will injure the cells which are responsible for hair growth by way of thermolysis (heat in the follicle tissue). This trauma is just below the threshold of injuring other cells (notably skin) around the follicle. IPL works best on clients who have light skin combined with dark hair growth. Should a hair follicle survive the IPL treatment, it will grow back slower and thinner. After a series of treatments the continual 'site specific' localized trauma to the follicle tissue will eventually destroy the hair growth. Once destroyed, a follicle cannot produce another hair.

Discomfort and Client Reaction


Remarkably, the sensation felt with an IPL treatment is no more intense than that of a laser (although much more total energy is being delivered to the tissue). The first two pulses from the IPL generally are not felt by the client. The third (if the technician chooses to deliver that many) begins to feel like a snap of a rubber band (as the tissue absorbs and retains the energy). Five or more pulses in less than one second may be uncomfortable for some clients. The use of a topical desensitizing spray or pre-treatment cold pack will alleviate most sensation.

Pre-Treatment Preparation

First spray the area with a lidocaine or benzocaine topical numbing compound. Now place a cold pack over the area intended to be treated for a minimum of 5 minutes prior to treatment. The patient will be able to withstand as many as 7 pulses per 2.5 cm²/second without much sensation. It is generally not necessary to deliver more than 5 pulses for effective hair removal treatment.

IPL Hair Removal Variation 1


This procedure is the most popular among professionals as it requires only shaving prior to the application. It will be necessary to remove all hair above the skin with a razor. Any hair which protrudes above the skin may cause pits in your optics which will reduce the overall efficiency of the device.



Turn your IPL driver to setting 6 pulse (54 jcm² per second) a place the IPL condenser against the skin in the area where the unwanted hair growth is occurring. Press the trigger switch then release (delivering four total emission equating to 54 jcm² of energy density) then move the condenser to the next treatment area. The treatment will effectively damage or destroy all anagen follicles in an area of 2.5 cm². An area 25cm square (about 10 inches) will required 100 pulses spaced evenly for effective hair removal.

IPL Hair Removal Variation 2: Deep Tissue Trauma With Follicle Extraction with Carbon Dye

Before applying treatment, remove all hair from the area by tweezing or waxing. IPL hair removal is most effective when applied to an empty follicle shaft. Human hair simply does not normally have enough pigment to allow for sufficient heat exchange to cauterize, desiccate and necrotize the cells which produce hair. To compensate for this lack of 'quantitative' and 'qualitative' photon targets, it will be necessary to place a high-density carbon dye inside the follicle prior to treatment.



Using a cotton-tipped applicator, completely cover the treatment area with the special dye included in your kit. Massage the dye into the follicle pore with a firm downward circular motion. Repeat 2-3 times to saturate the follicle pore. Use an **ethyl alcohol** based wipe (Isopropyl alcohol will not dissolve the dye) to lightly clean the excess from the surface of the skin. At this point you will have all desired follicles **visibly highlighted** with a dark spot (as seen above). Perform the standard IPL treatment as outlined in 'IPL Hair Removal Variation 1' as normal. *Don't over treat.* Applying 10 or more pulses in 1 second to the same area may result in a severe burn (leading to blistering).